

Vista



# Vista D

## Because life's best conversations happen on mobile phones too

### Made for all phones

Vista D R Li, Vista D R 312 and Vista D P Li hearing instruments connect directly to all phones – not just iPhones®

### Hands-free

There's no need to touch the phone. The patient clearly hears what the caller is saying, while microphones transmit their voice to their mobile phone.

### No streamer

Vista D R Li, Vista D R 312 and Vista D P Li acts like a Bluetooth® headset, connecting to mobile phones up to 10m away. Patients experience amazing sound quality and more enjoyable conversations.

### Easy to use

Patient needs were the top priority when designing Vista D R Li, Vista D R 312 and Vista D P Li, so it's easy to get them back into all their conversations, whether on a mobile phone or in person.



## Vista D direct connectivity



Vista D R Li   Vista D R 312   Vista D P Li

**Vista D R 312** Connects directly to all mobile phones\* – not just iPhones. **Vista D R Li** Direct connectivity in a versatile rechargeable hearing instrument. **Vista D P Li** First BTE to combine direct connectivity with lithium-ion rechargeable technology

# Your perfect style. Choose from products available on the new Vista D platform.

## Vista D P Li | **New**

### Behind-the-ear (BTE)

- › Stylish design
- › Easy to use Li-ion rechargeable
- › Made for all phones\*
- › Hands-free phone\* and video calls
- › Phone and stereo media audio streamed to both ears



## Vista D R Li | **with new telecoil option**

### Receiver-in-canal (RIC)

- › Stylish design
- › Telecoil option
- › Easy to use Li-ion rechargeable
- › Made for all phones\*
- › Hands-free phone\* and video calls
- › Phone and stereo media audio streamed to both ears



## Vista D R 312

### Receiver-in-canal (RIC)

- › Small, discreet design
- › Made for all phones\*
- › Hands-free phone\* and video calls
- › Phone and stereo media audio streamed to both ears



\*Smartphones and traditional mobile phones with a compatible Bluetooth Hands-Free Profile.

## Accessories

### TV Connector

Makes it simple to hear your favorite programs by wirelessly connecting your Vista D hearing aids to your TV. Just plug it in, press one button and you're all set to enjoy.



### Charger

Quickly charges your rechargeable hearing aids.



BTE charger

RIC charger

### Power pack

Charges your Vista D RIC rechargeable hearing aids when you are on the move.



### Hearing Remote app

Use your Apple or Android smartphone as a remote control for your Vista D hearing aids.



# Vista D RIC and BTE

		D 950	D 750	D 550	D 350	
<b>Soundscape Manager</b>						
<b>Auto Sound Control 3.0</b> Automatic program identifies and classifies sounds across environments, with a focus on conversations.	Music	•	•			
	Noise	•	•	•		
	Conversation in a crowd	•				
	Conversation in a small group	•	•			
	Conversation in quiet	•	•	•	•	
	Conversation in noise	•	•	•	•	
	Quiet	•	•	•		
	Media Control speech	•	•	•	•	
	Media Control music	•	•	•	•	
	Total streaming environments	2	2	2	2	
	Total listening environments	7	6	4	2	
	<b>Speech Target Pro</b> Provides optimal speech understanding and localization, even in the most challenging environments.	Speech Target Pro	•			
		Speech Finder	•	•		
		Speech Lock	•			
Speech Mapping		•				
Speech Target 2			•			
<b>Sound Director</b> Dynamically balances features to provide speech understanding, comfort and natural sound.	Speech intensifier	•	•	•	•	
	Noise reduction	•	•	•	•	
	Adaptive directionality	•	•	•	•	
	Sound Mapping	•	•	•		
<b>Sound Mapping</b> Makes it easier to localize sound, providing a more natural, realistic listening experience.	Personalized	•				
	Pinna Effect	•	•	•	•	
<b>Sound Suite</b>						
	Pulse Protector 2	•	•	•	•	
	Wind control	•	•	•	•	
	Feedback manager	•	•	•	•	
	Natural sound	•	•	•	•	
<b>Experience innovations</b>						
	Vista: trial	•	•	•	•	
	Vista: upgrade		•	•	•	
	Capture All	•	•	•	•	
<b>Convenience</b>						
	Direct connectivity	•	•	•	•	
	Rechargeable option	•	•	•	•	
	Wireless Synchronization	•	•	•	•	
	TV Connector	•	•	•	•	
	Binaural Phone	•	•	•	•	
<b>Fitting</b>						
	Automatic Adaptation Manager	•	•	•	•	
	Music Equalizer	•	•	•	•	
	Frequency compression 2	•	•	•	•	
	Tinnitus masker	•	•	•	•	
	IntelliVent	•	•	•	•	
	New first fit approach	•	•	•	•	
	Manual programs	•	•	•	•	
	Fitting channels	20	16	12	8	

## Fitting Range

