

Vista DX

Vista DX Product Portfolio 2.0



Vista

Vista DX

Superior speech understanding




Vista DX technology allows clients to have more speech understanding than ever before. Our new sound features help a client's ability to perceive where speech is coming from as well as the meaning behind words*. And now the Vista DX platform has an expanded product lineup featuring three new ITE styles and the new Vista DX M 312 BTE.

Enjoy increased flexibility in how you service clients with **remote adjust** – a new option that helps you deliver fine-tuning adjustments remotely.

*Singh, Kreuger, Besser, Wietoska, Launer, Meis. (2018). A Pre-Post Intervention Study of Hearing Aid Amplification: Results of the Emotional Communication in Hearing Questionnaire (EMO-CHeQ), IHCON 2018.

The intelligent sound processing of Vista DX helps enhance the subtle nuances of speech for improved understanding




Soundscape Manager provides more than the ability to hear sounds. It delivers details of speech to allow a more holistic hearing experience. The powerful, integrated Soundscape Manager makes Vista DX one of the world's most intelligent platforms, and it helps your clients to interpret more than just the speech they are hearing.

Identifies environments as accurately as a 20-year-old

Clients can hear their best across all listening situations with **Auto Sound control 3.0** and can enjoy life-enriching media with the help of Media Control. Years of training our algorithms with machine learning means this feature can identify life's ever-changing soundscapes just like a young, normal-hearing listener would.*

A well-balanced system

Provide clear speech, reduced noise, and overall sound quality that's comfortably balanced in any environment with **Sound Director**. And, our new **soft speech intensifier** feature emphasizes the details of speech to make it easier for clients to understand soft speech in challenging environments.



Hear speech with amazing clarity – even in noise

Address your clients' number one need to hear better in noisy environments with **Speech Target Pro**.** It uses binaural cues to know where speech is coming from and enhances it beyond what a typical hearing instrument can do. A new fourth dimension of Speech Target Pro, **Speech Separator**, was designed to help clients better understand words and the true meaning behind them.

Natural localization

Clients enjoy a realistic sound experience because sounds can be acoustically separated and perceived in their accurate locations. This is thanks to **Sound Mapping**, which brings back natural localization cues with an innovative 4 microphone strategy.

*Dr. David A. Eddins, WCA October 2018 Cape Town, South Africa.

**Ozmeral, Hoover, Gabbidon, and Eddins (2018). Development of Continuous Number Identification Test (CNIT). Ear and Hearing submitted for publication.

Incredible
begins with
the first words
they hear



And getting them
started has never
been easier

**A great first impression
in 5 simple steps**

Our Aura:fit fitting software enables a quick, easy and accurate fitting experience in 5 clicks or less. It ensures all the necessary settings are optimized, giving you confidence you're providing the best possible first fit with a better sounding starting point.

**Sound that eases them in
and keeps getting better**

Starting at a full target fit can be loud, harsh, and overwhelming, which can cause clients to stop wearing their hearing instruments. That's where our Automatic Adaption Manager (AAM) comes in. AAM starts with a volume and transparency that's comfortable to new wearers, then automatically transitions to peak performance so they gradually and seamlessly adjust to their new sounds.

More options to suit your clients

An expanded product lineup based on award-winning designs, and the ability to perform professional fine-tuning adjustments remotely.

Made for all connectivity

Your clients can enjoy easy direct connections to the people and content they love with hands-free phone calls, video calls and stereo media streaming to both ears – on all phones, with no intermediary accessories needed.*

SWORD 3.0 makes it possible

The Vista DX platform leverages the SWORD™ 3.0 (Sonova Wireless One Radio Digital) chip to provide made for all (MFA) connectivity including truly hands-free phone calls and media to both ears.


Get connected to 4 devices

Vista DX offers an extra level of convenience. Now, clients can use 2 Bluetooth® devices (phones, tablets) as well as 2 TV Connectors with their hearing instruments.



*Smartphones and traditional mobile phones with a compatible Bluetooth Hands-Free Profile. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

Remote adjust



Remote adjust gives you enhanced flexibility with the option to perform fine-tuning adjustments remotely. It offers you the ability to deliver professional care that's convenient for you and your clients, whether they're at home or on the move. Adjustments are easy to make and sent directly to your clients via the Hearing Remote app:

Clients can now enjoy professional care without having to go into your clinic. In just a few taps on the Hearing Remote app, they can apply your fine-tuning adjustments and hear the improved sound.



Give your clients more control

With the Hearing Remote* app your clients can also:


Change programs

Adjust volume

Mute and unmute

Adjust the balance between media and the environment

With the equalizer, your clients will have some level of personalization allowing them to adjust bass, mid and treble frequencies in the moment. All of these adjustments are captured in our Aura:fit software, so you have a complete picture of their preferences.



*For use of the Hearing Remote app, Vista hearing aids with Bluetooth connectivity are required. The Hearing Remote app works on Apple smartphones with iOS 12 or newer and Android smartphones with version 7 or newer. The equalizer screen may not appear as exactly as shown.

Fantastic hearing runs in the family

Say hello to our new family of stylish, comfortable hearing instruments, including an ITE with made for all (MFA) connectivity, and lithium m-ion rechargeable options proven to last all day. Plus, they're quick to charge and turn on and off automatically when removed from or inserted into the charger.

RIC and BTE



Vista DX R Li

Our smallest lithium-ion rechargeable Receiver-in-canal

Vista DX R Li T

Rechargeable Receiver-in-canal with telecoil

Vista DX R 312

Traditional battery Receiver-in-canal

NEW Vista DX M 312

Traditional battery Behind-the-ear

Vista DX P Li

Rechargeable Behind-the-ear

RIC and BTE hearing instrument colors

- Amber suede (S2)
- Espresso (P4)
- Platinum (P6)
- Beige (o1)
- Amber (P2)
- Pewter (P7)



Vista DX RIC charger



Vista DX P Li charger

ITE



NEW Vista DX W 312 Directional

Our first in-the-ear with MFA connectivity. Traditional battery.



NEW Vista DX 312 Omni

In-the-ear. Telecoil. Traditional battery.



NEW Vista DX 10A Omni


Our smallest in-the-ear. Telecoil. Traditional battery.

ITE hearing instrument colors


- Beige (26)
- Cocoa (22)
- Tan (14)
- Brown (28)
- Black* (o6)

Vista DX fitting ranges


RIC



BTE



ITE



Solutions for every need



PartnerMic

Lapel worn microphone for one-to-one conversations in noisy environments – streaming speech directly to your clients' hearing instruments.



Remote Control

Lets your clients discreetly adjust volume and switch programs at their convenience.

TV Connector

Clients can watch their favorite shows and movies in high-quality stereo sound by wirelessly connecting their hearing instruments to the TV or other media devices.



Roger Pen™



Roger Select™



Roger Table Mic II



RogerDirect™

Clients can benefit from the proven performance of Roger™ in noise and over distance; with RogerDirect, which allows Roger microphones to stream directly to Vista DX hearing instruments.




Products are not shown at actual size.

*Thibodeau, L. (2014). Comparison of speech recognition with adaptive digital and FM wireless technology by listeners who use hearing aids. American Journal of Audiology, 23(2), 201-210.

Combine it all with Vista:trial and Vista:upgrade for an unrivaled solution

Our Vista DX platform, along with Vista:trial and Vista:upgrade, can transform the hearing journey for you and your clients. Imagine overcoming typical first fit barriers by allowing clients to try hearing instruments for free in their daily lives. You can offer this with no commitment, meaning less pressure for them, and less stress for you. Give your clients the convenience they desire with the natural sound and speech clarity they expect.





Vista DX

Incredible sound performance to improve understanding of meaning in conversations

Proven lithium-ion rechargeable technology that lasts all day

Quick, easy and accurate first fits with sound that eases them in and keeps getting better

Enhanced flexibility with remote fine-tuning adjustments

A sleek and comfortable full portfolio of products, including an ITE with made for all (MFA) connectivity

Hands-free phone calls, video calls and stereo media streaming to both ears, on all phones

Risk-free evaluation in their everyday lives with Vista:trial

Future-proof hearing instruments with Vista:upgrade

At a glance

Soundscape Manager		Vista DX 950	Vista DX 750	Vista DX 550	Vista DX 350
Environmental classification	Auto Sound control 3.0	•	•	•	•
	Total listening environments	7	6	4	2
	Conversation in a crowd	•			
	Conversation in a small group	•	•		
	Music	•	•		
	Quiet	•	•	•	
	Noise	•	•	•	
	Conversation in quiet	•	•	•	•
	Conversation in noise	•	•	•	•
	Total streaming environments	2	2	2	2
	Media Control Speech	•	•	•	•
	Media Control Music	•	•	•	•
Sound optimization	Sound Director	•	•	•	•
	Speech intensifier	•	•	•	•
	Soft speech intensifier	•	•	•	
	Noise reduction	•	•	•	•
	Microphone features	•	•	•	•
Localization	Sound Mapping	•	•	•	
	Personalized	•			
	Pinna Effect	•	•	•	•
Performance in challenging environments	Speech Target Pro	•			
	Speech Finder	•	•		
	Speech Lock	•			
	Speech Mapping	•			
	Speech Separator	•			
	Speech Target 2		•		
Fine-tuning channels		20	16	14	12

Available in all technology levels

Sound stabilization

- › Pulse protector 2
- › Wind control
- › Feedback manager
- › Natural sound

Fitting

- › New first fit approach
- › Automatic Adaptation Manager
- › Frequency compression 2
- › Tinnitus masker
- › Music Equalizer
- › Manual programs
- › IntelliVent

Ease and convenience

- › Remote adjust
- › Made for all direct connectivity
- › TV Connector
- › Remote Control
- › PartnerMic
- › Hearing Remote app
 - Equalizer
- › Wireless synchronization
- › Binaural Phone
- › RogerDirect
 - Roger Pen
 - Roger Select
 - Roger Table Mic II

Experience innovations

- › Vista:trial
- › Vista:upgrade
- › Capture All

Compatibility may be limited by product style